



Bio-Therapeutic
Skin Care for life inside and out®

bt-zoom® quick start guide



1. Read manual prior to using the equipment.
2. Open battery compartment on base of handle by gently pressing **OPEN** latch.
3. Install 3 AAA batteries. To replace battery cover, insert small tab at top of battery compartment and push in.

4. Press the **POWER BUTTON** once to turn on lights. Continue to press power button to access light options and to turn off.





5. If using silicone sleeve, place handle of bt-zoom® in corresponding area of sleeve and stretch around lens.